Cardinal Health™ technical data sheet

Kendall™ ECG Electrodes H92SG

Ref. Code: 31.1925.21

Product information

Shape / size	Oval / 57 x 34 mm
Total product surface (incl. grip)	1522 mm ²
Gel area	254 mm ²
Adhesive area	1050 mm ²
Sensor area	50 mm ²
Product thickness (adapter excluded)	1 mm
Adapter	Stud
Integrated lead wire (length / colour)	-

Application

Long term Monitoring

Materials information

Backing material	Polyethylene foam (PE), white	
Adhesive characteristics	Medical grade acrylic-based pressure sensitive adhesive	
Gel characteristics	Conductive and adhesive acrylic-based hydrogel	
Supporting / back label	Polypropylene film (PP), transparent	
Release liner	PE/Paper/PE, one side siliconized, white	
Sensor	Polymer Ag/AgCl coated	
Adapter / connector	Stainless steel	
Integrated lead wire jacketing	-	
Integrated lead wire cord	-	
•••••••••••••••••••••••••••••••••••••••		

Electrical Characteristics

(ANSI/AAMI EC 12, average measured before packaging)

ACZ I impedance (before defibrillation simulation)	200 Ohm
DC offset voltage (before defibrillation simulation)	0.2 mV
SDR (remaining potential after defibrillation simulation)	11 mV
Slope (potential decline after defibrillation simulation)	0.2 mV/s
COIIN (combined offset instability and inner noise)	4 μV
Bias current tolerance (DC offset voltage after DC loading)	6 mV





Cardinal Health™ technical data sheet Kendall™ ECG Electrodes H92SG

Special Features

MRI compatibility (≤ 3T)	No
X-ray translucence	No
Peel off resistance, max. (from PE foil, 180 °, 50 mm/min, ASTM D1876)	5 N
Integrated abrader	No
Repositionability	No

Biocompatibility

Test acc. to DIN ISO 10993	Passed
LATEX content	Latex free
PVC content	PVC free
DEHP content	DEHP free

Environment

RoHS directive	In compliance	

Packaging

	L x W x H (mm)	Weight (Kg)	Quantity (electrodes)	Material
Product packaging (pouch)	220 x 180	-	50	Paper / PE / Alu / PE
Department packaging (box)	194 x 194 x 240	0.72	500	Paper cardboard
Box recycled	No			

	-	
Box recycable	Yes	
Box bleached	Yes	

Shelf life

Keep dry

Product shelf life in an intact, originally sealed packaging

Do not use if package is opened or damaged

24 months between 10-30°C

CE

Keep away from sunlight

Legal Status

CE Mark

Class I Non-Sterile

cardinalhealth.co.uk

For healthcare professionals only.

Important information: Before using any medical device, review all relevant package information, including the label and/or the Instructions for Use. CARDINAL HEALTH, Cardinal Health LOGO, ESSENTIAL TO CARE and KENDALL are trademarks of Cardinal Health and may be registered in the US and/or in other countries. © 2021 Cardinal Health. All Rights Reserved. Lit. No. 2GM19-985745-2 12/2021

